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C006-3 F1:5ml×2
F2:5ml×2

Explanation of Symbols

	Sterilized using aseptic processing techniques
	Batch code
	Date of manufacture
	Use-by date
	Temperature limitation
	Keep away from sunlight

Precautions and Warnings

- Before using this product, users should read and understand the operation and precautions.
- This product cannot be used for injection.
- This product is sterilized by filtration and is not suitable for re-sterilization.
- Do not use this product if:
 - Product packaging is damaged, or the seal has been destroyed.
 - The product exceeds the shelf life.
 - The product has become discolored, cloudy, or granular.
- This product contains gentamicin sulfate. Appropriate precautions should be taken to ensure that the patients are not allergic to the antibiotic.
- This product contains human serum albumin derived from human blood. The human serum albumin in the product has been tested to be negative for HIV-AB, TP and anti-HCV, and HBsAg, HBsAb, HBeAg, HBeAb and HBcAb all meet the requirements. All blood products should be considered potentially infectious.
- This product is suitable for single use. After opening the bottle, discard the remaining unused product.

8. In order to avoid contamination, please use aseptic technique.

9. This product can be used with Kitazato products of mainstream companies in the market, mouse embryo assay(MEA)-1-Cell MEA≥80%, the test results are as follows:

Brand	Vitrification thawing lot number	Recovery rate	Blastocyst formation rate
Kitazato	Y20060201	100%	94.4%

Other manufacturers' s products have not been tested.

10. The product shall remain stable and effective within the specified storage conditions and the shelf life.

11. This product contains dimethyl sulfoxide, which is detected according to the national standard GB/T 16886.5-2017 biological Evaluation of Medical Devices part 5: In Vitro Cytotoxicity Test.

Intended use

This product is suitable for vitrification and freezing of embryo.

Composition

L - serine	Taurine
L - valine	Glycine
L - leucine	Glucose
L - proline	Sodium lactate
L - threonine	Ethylene glycol
L - isoleucine	Calcium chloride
L - methionine	Sodium chloride
L - tryptophan	Magnesium sulfate
L - glutamic acid	Potassium chloride
Sodium hydroxide	
L - arginine hydrochloride	
L - cystine dihydrochloride	
L - tyrosine disodium salt dihydrate	
4 - hydroxyethyl piperazine ethyl sulfonic acid	

L - asparagine
L - aspartic acid
L - phenylalanine
Dimethyl sulfoxide
Pyruvic acid sodium
Sodium bicarbonate
N(2) - L - alanyl - L - glutamine
L - lysine monohydrochloride
L - histidine hydrochloride hydrate
Potassium dihydrogen phosphate
Human Serum Albumin(10mg/mL)
3 - (N-morpholine) propyl sulfonic acid

F1

Base ingredients: Gentamicin sulfate (7.5µg/mL).

F2

Basic ingredients: Gentamicin sulfate (5.4µg/mL), trehalose.

Storage and Transportation

- Store dark at 2°C~8°C for 6 months.
- Store dark at 2°C~8°C during transportation.

Product performance index

- 1.pH (37°C) : 7.2~7.6 .
2.Osmolality :
F1: 222 mOsmol/kg ~ 422 mOsmol/kg (diluted 10 times);
F2 : 494 mOsmol/kg ~ 694 mOsmol/kg (diluted 10 times).
3. Endotoxin : < 0.2EU /mL.
4. Product sterility : Pass.
5. 1-Cell Mouse embryo assay : ≥80%.

Contraindications

Those who are allergic to human serum albumin are prohibited.

Operation

- Procedures are to be done at room temperature (22-28°C) for at least 30 minutes to use.
 - Prepare a dish with 1.0mL of F1.
 - Carefully transfer the embryos (no more than three at a time) to balance in F1 for 6 minutes to a maximum of 15 minutes, with minimal removal of culture medium.
 - When the embryos were balanced in F1, 1.0mL F2 was taken and added to the dish.
 - Carefully remove the embryos equilibrated in F1 to F2, under microscope observation, gently pipette the embryos, if floating, and place at the bottom. With minimal volume less than 1µL of F2, and immediately plunge it into LN(2). Whole procedure should be completed in 50-70s.
- Note:
- If vitrification is required to freeze more embryos, repeat steps 2-5 above, using new utensils and fresh frozen liquid 1(F1) and frozen liquid 2(F2).
 - When freezing blastocyst, artificial shrinkage should be performed first, and then according to steps 1~5. Try to freeze the blastocyst of grade 4 or above, which is the stage of expansion.
 - Try to avoid bubbles when handling all liquids.
 - Minimize the carry-over of the culture medium each time the embryo is transferred.
 - F1 and F2 should be operated at room temperature.