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C007-4 W1:5ml×2 W2:5ml×1  
W3:5ml×1 W4:5ml×1

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C007-4 W1:5ml×2 W2:5ml×1  
W3:5ml×1 W4:5ml×1

### Explanation of Symbols

	Sterilized using aseptic processing techniques
	Batch code
	Date of manufacture
	Use-by date
	Temperature limitation
	Keep away from sunlight

### Precautions and Warnings

1. Before using this product, users should read and understand the operation and precautions.
2. This product cannot be used for injection.
3. This product is sterilized by filtration and is not suitable for re-sterilization.
4. Do not use this product if:
  - Product packaging is damaged, or the seal has been destroyed.
  - The product exceeds the shelf life.
  - The product has become discolored, cloudy, or granular.
5. This product contains gentamicin sulfate. Appropriate precautions should be taken to ensure that the patients are not allergic to the antibiotic.
6. This product contains human serum albumin derived from human blood. The human serum albumin in the product has been tested to be negative for HIV-AB, TP and anti-HCV, and HBsAg, HBsAb, HBeAg, HBeAb and HBcAb all meet the requirements. All blood products should be considered potentially infectious.
7. This product is suitable for single use. After opening the bottle, discard the remaining unused product.

8. In order to avoid contamination, please use aseptic technique.

9. This product can be used with Kitazato products of mainstream companies in the market, mouse embryo assay(MEA)-1-Cell MEA $\geq$ 80%, the test results are as follows:

Brand	Vitrification thawing lot number	Recovery rate	Blastocyst formation rate
Kitazato	Y20060901	100%	96.3%

Other manufacturers' products have not been tested.

10. The product remain stable and effective within the specified storage conditions and the shelf life.

### Intended use

This product is suitable for thawing and rewarming of vitrified embryos.

### Composition

L - valine  
L - serine  
L - proline  
L - leucine  
L - threonine  
L - isoleucine  
L - asparagine  
L - tryptophan  
L - methionine  
L - aspartic acid  
L - glutamic acid  
L - phenylalanine  
Pyruvic acid sodium  
L - cystine dihydrochloride  
N(2) - L-alanyl -L-glutamine  
L - lysine monohydrochloride  
L - arginine monohydrochloride  
3 - (N-morpholine) propionic acid  
Potassium dihydrogen phosphate  
L - histidine hydrochloride hydrate  
Human Ssrum Albumin(10mg/mL)  
L - tyrosine disodium salt dihydrate  
4-hydroxyethyl piperazine ethyl sulfonic acid

W1:

Base ingredients: Trehalose, Gentamicin sulfate (7.1 $\mu$ g/mL).

W2:

Base ingredients: Trehalose, Gentamicin sulfate (8.1 $\mu$ g/mL).

W3:

Base ingredients: Trehalose, Gentamicin sulfate (8.7 $\mu$ g/mL).

W4:

Base ingredients: Gentamicin sulfate (9.2 $\mu$ g/mL).

### Storage and Transportation

1. Store dark at 2°C~8°C for 6 months.
2. Store dark at 2°C~8°C during transportation.

### Product performance index

1.pH (37°C) : 7.2~7.6 .

2.Osmolality :

W1: 109 mOsmol/kg~169 mOsmol/kg (diluted 10 times);

W2 : 63 mOsmol/kg to 94 mOsmol/kg (diluted 10 times);

W3 : 521mOsmol/kg~621mOsmol/kg ;

W4 : 255mOsmol/kg~280mOsmol/kg .

3. Endotoxin : < 0.2EU /mL.

4. Sterility testing : Pass.

5. 1-Cell Mouse embryo assay :  $\geq$ 80%.

### Contraindications

Those who are allergic to human serum albumin are prohibited.

### Operation

1. Bring the quantity to be used of W1 and equilibrate for at least 40 minutes at 37°C, W2, W3 and W4 to room temperature (22-28°C) for at least 30minutes prior to warming vitrified specimens.

2. Dispense a minimum volume of 1mL of W1 and put it on 37°C heated stage.

3. Prepare a dish with 500 $\mu$ L droplets of W2, W3, W4 and place the dish at room temperature.

4. Remove the frozen carrier which with the embryo facing up, and immediately plunge it into W1. Under microscopic observation, gently move the carrier until the embryos are released, and gently pipette the embryos, if floating, and place at the bottom. Leave the embryos for a total of 50-60s in the W1.(This operation needs to be carried out at 37°C)

5. Transfer the embryos into W2 and remain for 3 minutes (This operation needs to be carried out at room temperature).

6. Transfer the embryos into W3 and remain for 5 minutes (This operation needs to be carried out at room temperature).

7. Transfer the embryos into W4 and remain for 1 minute, then transfer it to culture medium.

Operation such as assisted incubation can be completed in W4. (This operation needs to be carried out at room temperature)

Note:

① If more embryos need to be rewarmed and thawing, repeat the above steps 2 to 7, using new utensils and fresh thawing solution.

② When operating all culture medium, avoid bubbles as far as possible.

③ Every time you transfer an embryo, you should try to reduce the amount of culture fluid.

④ W1 is operated at 37°C; W2 to are operated at room temperature.