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住所：四川省成都市天府新区新兴街道天工大道916号A2栋3-4层  
Residence: 3rd-4th floor, Building A2, Tiangong Road No. 916, Xinxing street, Chengdu, Sichuan Province People's Republic of China

生产地址：四川省成都市天府新区新兴街道天工大道916号A2栋3-4层  
Manufacturing address: 3rd-4th floor, Building A2, Tiangong Road No. 916, Xinxing street, Chengdu, Sichuan Province People's Republic of China

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## Explanation of Symbols

|  |  |
|--|--|
|  | Sterilized using aseptic processing techniques |
|  | Batch code                                     |
|  | Date of manufacture                            |
|  | Use-by date                                    |
|  | Temperature limitation                         |
|  | Keep away from sunlight                        |

## Precautions and Warnings

1. Before using this product, users should read and understand the Operation and precautions.
2. This product cannot be used for injection.
3. This product is sterilized by filtration and is not suitable for re-sterilization.
4. Do not use this product if:
  - The product package is damaged or the seal has been destroyed.
  - The product exceeds the shelf life
  - The product has become turbid or granular.
5. This product contains gentamicin sulfate. Appropriate preventive measures should be taken to ensure that the patients are not allergic to the antibiotic.
6. This product contains human serum albumin derived from human blood. The human serum albumin in this product has been tested to be negative for HIV-Ab, TP, and Anti-HCV, and HBsAg, HBsAb, HBeAg, HBeAb, HBCAb all meet the requirements. All blood products should be considered

potentially infectious. Appropriate precautions should be taken to ensure that the patients do not have an allergic reaction to the human serum albumin.

7. In order to avoid contamination, the media should be used only with aseptic technique.
8. After taken out, the product cannot be poured back to the vial for further storage and use.
9. The product remains stable and effective within the specified storage conditions and the shelf life.
10. Caution: Federal law restricts this device to sale by or on the order of a physician.

## Intended use

This product is intended for cryopreservation of human semen.

## Composition

|                                |                    |
|--------------------------------|--------------------|
| L-arginine                     | L-Serine           |
| L-cystine                      | Sodium chloride    |
| L-phenylalanine                | Sodium pyruvate    |
| L-asparagine                   | Potassium chloride |
| L-histidine                    | Calcium chloride   |
| L-aspartic acid                | HEPES              |
| L-isoleucine                   | MOPS               |
| L-glutamic acid                | Sodium hydroxide   |
| L-leucine                      | Glycerol           |
| L-methionine                   | Taurine            |
| L-proline                      | Magnesium sulfate  |
| L-threonine                    | Sodium bicarbonate |
| L-tryptophan                   | Glycine            |
| L-tyrosine                     | Glucose            |
| L-valine                       | Sodium lactate     |
| L-lysine monohydrochloride     |                    |
| N (2) -L-alanyl-L-glutamine    |                    |
| Human Serum Albumin (5mg / mL) |                    |
| Potassium dihydrogen phosphate |                    |
| Gentamycin sulfate (10 µg/mL)  |                    |

## Storage and Transportation

1. Store dark at 2°C ~ 8°C for 6 months.
2. Store dark at 2°C ~ 8°C during transportation.
3. Product can be used within 7 days from first opening when store dark at 2°C ~ 8°C.
4. The remaining media should be discarded 7 days after opening the bottle.

## Product performance index

1. pH(37°C): 7.2 ~ 7.6.
2. Osmolality:  
200 mOsmol / kg ~ 280 mOsmol / kg (diluted 10 times)
3. Product sterility: Pass.

## Contraindications

Those who are allergic to human serum albumin and gentamicin sulfate are prohibited.

## Operation

1. Product preparation  
Pre-warm the corresponding volume of Sperm Freezing Medium at room temperature (20-25°C) to ensure that the media reaches temperature equilibrium (at least 30 minutes).
2. Semen preparation  
The semen is allowed to liquefy for 37°C at least 30 minutes.
3. Freeze semen  
Transfer the liquefied semen to a sterile 15 mL conical centrifuge tube. Dilute the liquefied semen with equal volume of Sperm Freezing Medium in 30-40s. Add the medium drop by drop onto the semen and mixed carefully after each addition. Aliquot the semen mixture into sterile and non-enzymatic cryotubes.
- 3.2 Keep cryotubes at room temperature for about 10 minutes before freezing.
- 3.3 Transfer the cryotubes quickly into the liquid nitrogen and store at -196 °C.  
\* Suspend cryotubes about 5 cm above the surface of the liquid nitrogen for 10 minutes. Transfer the cryotubes into the liquid nitrogen

and store.  
2. Thaw semen  
Remove cryotubes from -196 °C and place them in a water bath at 37 °C for 5-10 minutes.

## Sperm Freezing Medium